Trimercaptotriazine (TMT)

Supra Sciences

Trimercaptotriazine Resin, MP, is a macroporous polystyrene resin functionalized with a trimercaptotriazine end group while its Silica counterpart is bound to a silica backbone. This reagent has shown to be an exceptional scavenger of organometallic reagents. MP-TMT is also extremely efficient in its ability to scavenge Copper.

S N SH

The limitations of solution phase trimercaptotriazine are based primarily around solubility in various solvents. The macroporous backbone of MP-TMT allows for its use in all solvents, and the Silica counterpart is effective in all solvents except any protic media. Simple filtration removes it (and the complexed metal) from your reaction.

General Reaction

References

Rosso, V.W. *Org. Process Res & Dev.* **1997**, *1*, 311-314 Ishihara, K. *Chem Lett,* **2000**, *29*, 1218-1219 Chen,C;Dagneau, P. *J.Org. Chem.* **2003**, *68*, 2633-2668.

Solvent Compatibility

MP:	THF	Silica:	THF
	DMF		DMF
	NMP		NMP
	DCM		DCM
	DCE		DCE
	MeOH		ACN
	EtOH		

Ordering Information

MP-Trimercaptotriazine

Loading: 3.0-3.4 mmol/g	10g	SPMP 26-10
	25g	SPMP 26-25
Bead size: 330-1225 microns, 15-50 mesh	100g	SPMP 26-100
(>90% within)	1Kg	SPMP 26-1kg

Trimercaptotriazine-Fine

Loading: 1.7-1.8 mmol/g	10g	SPMP 27-10
	25g	SPMP 27-25
Bead size: 100-200 mesh	100g	SPMP 27-100
	1Kg	SPMP 27-1kg

Si-Trimercaptotriazine

Loading: 0.3-0.35 mmol/g	10g	SPSi 22-10
	25g	SPSi 22-25
Bead size: Avg. 40-62 microns	100g	SPSi 22-100
	1Kg	SPSi 22-1kg

For additional information contact info@suprasciences.com or visit www.suprasciences.com